

# geolometer

R O C K   S O L I D L Y   R E L I A B L E



**GL-AT 114**

**Low temperature constant  
temperature groove**

## PRODUCT APPLICATION

Low temperature constant temperature slot is one of the commonly used equipment in the laboratory, it is widely used in petroleum, chemical, metallurgy, medicine, biochemical testing, chemical analysis and research department, colleges and universities, factory laboratories and measurement quality inspection department, the equipment provides a hot and cold controlled, uniform temperature and constant maintenance environment for the sample.

## Product Specifications

SPECIFICATIONS	GL-AT 114
Temperature range	-6.00 to 59.99°C
Volatility	±0.1°C
Capriciousness	15L
Operating slot volume	3:250×200×300mm
Pump flow rate	3.6L/mm
Power	1.5kW
Power Supply	AC220V ± 10%, 50HZ